

STANDARD OPERATING PROCEDURES
LABORATORY

PROCEDURE TITLE: Preparing TFM and Niclosamide Solutions

SCOPE: To describe procedures for preparing various TFM and Niclosamide stock solutions and standards.

A. 99% pure TFM stock solution \approx 1000 mg/L

1. Record initial and final weight of the container, and weight of sample removed, on the appropriate form in the Chemical Register and User Log.
2. Tare a 100 mL volumetric flask.
3. Accurately weigh out 0.100 g of 99% TFM.
4. Record weight (X) to the nearest 0.1 mg.
5. Add approximately 3 mL Dimethylformamide (DMF) to dissolve TFM and dilute up to 100 mL using deionized or HPLC grade water.
6. Calculate concentration of TFM stock solution:

$$\text{mg/L stock TFM} = \frac{(X \text{ mg TFM})(.99)}{100 \text{ mL}} \left(\frac{1000 \text{ mL}}{1 \text{ L}} \right)$$

B. TFM working standards

1. Prepare working standards by diluting 99% pure TFM stock solution with test water to desired concentrations which bracket the expected concentration in the samples.
2. Use equation to determine X mL of stock solution required to prepare each concentration:

$$X \text{ mL stock TFM} = \frac{\text{desired Conc. TFM}(\text{mg / L}) * \text{desired Vol. Std.}(\text{mL})}{\text{conc. stock TFM}(\text{mg / L})}$$

C. Standards for TFM quality assurance testing

1. Prepare three TFM standards of about 700 mg/L, 1100 mg/L, and 1500 mg/L.
 - a. Tare 100 mL volumetric flask on the analytical balance.
 - b. Log 99% pure TFM weight data (initial container weight, sample weight, and final container weight) on the appropriate form in the Chemical Register and User Log.
 - c. In a tared 100 mL volumetric flask, weigh out about 0.07g, 0.11g, or 0.15g of 99% pure TFM, dissolve in approximately 3mL dimethylformamide (DMF) and dilute to 100.0 mL with HPLC grade water.
 - d. Record weights to the nearest 0.1 mg and calculate standard concentrations from the recorded weights.

$$mg/L\ TFM = \frac{(g\ TFM)(1000\ mg)(1000\ mL)}{(100\ mL)(g)(L)}$$

D. Niclosamide stock solution (100 mg/L)

1. Log container initial, sample, and container final weight data into the chemical log book, on the appropriate form.
2. In a 100 mL volumetric flask accurately weigh out 0.010 g of >99% niclosamide.
3. Record weight (X) to the nearest 0.1 mg.
4. Dissolve the niclosamide with DMF and dilute up to 100 mL using DMF.
5. Calculate concentration of niclosamide stock solution:

$$mg/L\ stock\ niclosamide = \frac{(X\ mg\ niclosamide)(.994)\left(\frac{1000\ mL}{1\ L}\right)}{100\ mL\ DMF}$$

E. Niclosamide field standard (100 mg/L) for Sea Lamprey Control Agents

1. In the 99% Niclosamide Log, be sure to log in the weight data (container initial, sample, and container final weight).
2. In a 50 mL beaker accurately weigh out 0.100 g of \geq 99% niclosamide.
3. Record weight (X) to the nearest 0.1 mg.

4. Dissolve the niclosamide in DMF, transfer to a 1000 mL volumetric flask, rinsing well with a wash bottle filled with DMF, and dilute up to 1000 mL using DMF.
5. Mix by inverting . Use the ultrasonic bath if necessary, to aid Niclosamide going into solution.

Approved by: _____ **Date:** _____
Laboratory Supervisor

Reviewed by: _____ **Date:** _____
Quality Assurance Representative